

# Certificate of Conformity

## Pharmatech Dust Cap Seal

Pharmatech confirms that the Dust Cap Seal enclosed has been manufactured from CFR compliant virgin, silicone rubber. This Dust Cap Seal has been inspected according to Pharmatech Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter: 150mm\*  
Batch No: 5039/7115  
Part No: 095/150/T  
Material: Silicone Rubber  
Colour: Translucent  
Curing: Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to CFR Standard:  
Title 21, Chapter 1, Schedule 177.2600  
U.S.P. Class VI

#### TSE/BSE Statement:

We confirm that no materials derived from Animal Origin are used in the manufacture of these products, and no intermediates and/or auxiliary agents of animal origin have been used.

# 150 T

Approved by: Rehinder Date: 12 MAY 16

Name: SHAWN O'LODAN Position: DOCUMENT CONTROLLER

**Pharmatech** (a division of Terriva Limited)

4-5 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@pharmatech.co.uk](mailto:info@pharmatech.co.uk)